

# AS 200 JET PRO & AS 200 JET PHARMA

## KEY ADVANTAGES



### | Innovative air jet sieving technology

Enables efficient and reproducible **particle size analysis** of fine powders, especially **below 100 µm**.

### | All-in-One Device → Sieving → Weighing → Evaluating

**Saves lab space** and **eliminates errors** from transferring material between devices. **Streamlines workflow** for efficiency.

### | Automatic underpressure regulation

Ensures consistent and reliable sieving results by **maintaining optimal air flow** throughout the process.

### | User-friendly integrated software and touch display

**Easy control of all parameters, fast evaluation, and reporting.** Improves user experience with **guided sieving** and process reliability.

### | Data export and connectivity

**USB** and **Ethernet interfaces** for seamless data transfer and integration into **Laboratory Information Management Systems**

### | Various checks and assistants

**Reduce operating error**, ensure reproducible results, and make operation easy even for new users.

### | Pharma-conform design

Meets the strict requirements of the pharmaceutical industry, including **full GMP compliance** and easy cleaning thanks to 316L steel sieving chamber. Includes **audit trail, user management** and **e-signatures**.

### | IQ/OQ documentation available

Supports qualification and validation processes, **simplifying regulatory compliance**.

*SMART.*

*PRECISE.*

*EFFICIENT.*

